

InSePtion®

Interspinous MIS Fixation Device

The InSePtion® MIS Posterior Spinal Fixation System combines truly unique implant functionality with ergonomic and intuitive surgical instrumentation. InSePtion offers a full range of implant options to accommodate varying patient anatomy.

InSePtion® is a registered trademark of the USPTO.




Kalitec Medical

www.kalitecmed.com


info@kalitecmed.com (407) 545-2063

14-MBR-0000 | REV_D 04/19

FEATURES




Up to 10° Transverse Angulation
Accommodates varying patient bony anatomy



Color Coded Implants

Sizes range from 8 to 18mm in 2mm increments.
All implants are manufactured from Titanium Alloy (Ti-6Al-4V ELI)



Kalitec Medical

www.kalitecmed.com

info@kalitecmed.com (407) 545-2063